



NUTRITIONIST PRO™ SYSTEM REQUIREMENTS

Thank you for purchasing Nutritionist Pro™. For your convenience, our system requirements for our product are as follows:

Single User Version Installation

- ✓ Windows 2000, XP, Vista & Windows 7
- ✓ Pentium processor 166 MHz (300 or better recommended)
- ✓ 64 MB RAM (128 MB RAM recommended)
- ✓ 200 MB available hard drive space to install
- ✓ CD-ROM drive
- ✓ SVGA monitor and video card that supports at least 800 x 600 pixels and 256 colors
- ✓ HP or compatible laser printer

Network Version: Client Workstation Installation

- ✓ Windows 2000, XP, Vista & Windows 7
- ✓ Pentium processor 166 MHz (300 or better recommended)
- ✓ 64 MB RAM (128 MB RAM recommended)
- ✓ 50 MB available hard drive space to install
- ✓ CD-ROM drive
- ✓ SVGA monitor and video card that supports at least 800 x 600 pixels and 256 colors
- ✓ HP or compatible laser printer

Network Version: Server Installation

- ✓ Windows Server 2000, 2003, 2008 or workstations acting as a server running XP Pro, Vista & Windows 7
- ✓ Our network version is compatible with VMware®
- ✓ Pentium processor 166 MHz (300 or better recommended)
- ✓ 128 MB RAM
- ✓ 350 MB available hard drive space to install
- ✓ CD-ROM drive
- ✓ SVGA monitor and video card that supports at least 800 x 600 pixels and 256 colors

Nutritionist Pro can be used on a MAC if a 3rd party emulator program is used.

For information about how to use Nutritionist Pro on a MAC, please visit www.nutritionistpro.com/technical_faq.php or call 1-800-709-2799 ext 9

